# Cochlear™ Vistafix®

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### INTERNATIONAL ORIGINAL ARTICLES

EUGENE N. MYERS, MD  
International Editor

#### Long-term study of quality and safety of osseointegration for the retention of auricular prostheses

THOMAS WESTIN, MD, PhD, ANDERS TJELLSTRÖM, MD, PhD, EVA HAMMERLID, MD, KERSTIN BERGSTRÖM, MD, and  
BO RANGERT, MD, Göteborg, Sweden

The aim of this article is to describe the safety and quality of the osseointegrated implant technique for the retention of craniofacial prostheses, to present a protocol for collection of clinical data, and to discuss the impact of the procedure on the patient quality of life. A protocol was designed and used to study patients who had received auricular prostheses consecutively since 1979 at our department. The patients were asked to answer a questionnaire designed to describe symptoms and problems specific for someone wearing an auricular prosthesis. In total, 99 patients received 107 prosthetic ears (8 patients had bilateral defects) retained on 309 implants (2 to 4 implants/ear). Patients of all ages were represented, and only 9 discontinuities were reported. Most patients (95%) wear their prosthesis every day, in most cases more than 10 hours/day. The follow-up period ranged from 1 to 12 years, giving a total of 2624 postoperative observations of implants, with a 3% incidence of significant skin reaction. We conclude that the surgical technique for auricular prostheses retained on osseointegrated implants is simple and associated with a low rate of perioperative and long-term complications. It offers a high degree of stability and aesthetic satisfaction. (*Otolaryngol Head Neck Surg* 1999;121:133-43.)

From the Department of Otolaryngology-Head and Neck Surgery, Sahlgrenska sjukhuset, University Hospital of Göteborg (Drs. Westin, Tjellström, Bergström, and Hammerlid), and Nobelpriset AB (Dr. Rangert). Presented at the Annual Meeting of the American Academy of Otolaryngology-Head and Neck Surgery, New Orleans, La., Sept. 17-20, 1995. Reprint requests: Associate Professor Thomas Westin, MD, PhD, Department of Otolaryngology-Head and Neck Surgery, Sahlgrenska sjukhuset, University of Göteborg, S-413 45 Göteborg, Sweden. Copyright © 1999 by the American Academy of Otolaryngology-Head and Neck Surgery Foundation, Inc. 0194-2998/99/95.00 + 0 23/7789627

Plastic surgical reconstruction of the defective ear is a major challenge.<sup>1</sup> In selected cases an auricular prosthesis presents an attractive and acceptable alternative. Modern techniques and plastic materials make it possible to make a prosthesis of such high quality that it is difficult to distinguish it from a normal ear.<sup>2</sup> Experience from the successful retention of dental prostheses on osseointegrated implants<sup>3</sup> has been used in developing a system for the retention of auricular prostheses. The surgical technique for the osseointegration method when used in craniofacial reconstructions is well known, and the results are well documented.<sup>4-6</sup> The biologic and physical properties of the implant system have been thoroughly investigated.<sup>7</sup>

The aim of this study was to design and present a protocol enabling systematic collection of clinical data during clinical follow-up visits to ensure the safety and quality of the system. This article describes the patients and the problems they have encountered after receiving auricular prostheses consecutively since 1979 at the Department of Otolaryngology-Head and Neck Surgery, Sahlgrenska Hospital, University of Göteborg, Sweden. The quality of the system was further investigated by a questionnaire designed to examine the way in which the prostheses affected the quality of life of the patients.

#### METHODS AND MATERIAL

##### Protocol


The protocol presented in this study was created to enable registration (and processing) of adequate data not only for the patients but also for each individual implant (or fixture). The data registered for the patients were considered to be of importance for the outcome of the prosthesis system. The patient history includes diagnosis, general history, and when relevant, tumor data. It also includes data on smoking, medication, and diseases, along with a section on prosthesis data. The protocol sheets are presented as an appendix.

The protocol shows a coordinate system in which each individual fixture inserted can be given a position. The coordinates are made clockwise, with the center in the ear canal.

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## INTERNATIONAL ORIGINAL ARTICLES

EUGENE N. MYERS, MD  
International Editor

### Long-term study of quality and safety of osseointegration for the retention of auricular prostheses

THOMAS WESTIN, MD, PhD, ANDERS TJELLSTRÖM, MD, PhD, EVA HAMMERLID, MD, KERSTIN BERGSTRÖM, MD, and  
BO RANGERT, PhD, Göteborg, Sweden

The aim of this article is to describe the safety and quality of the osseointegrated implant technique for the retention of craniofacial prostheses, to present a protocol for collection of clinical data, and to discuss the impact of the procedure on the patient quality of life. A protocol was designed and used to study patients who had received auricular prostheses consecutively since 1979 of our department. The patients were asked to answer a questionnaire designed to describe symptoms and problems specific for someone wearing an auricular prosthesis. In total, 99 patients received 107 prostheses (8 patients had bilateral defects) retained on 309 implants (2 to 4 implants/ear). Patients of all ages were represented, and only 9 discontinuities were reported. Most patients (95%) wear their prostheses every day, in most cases more than 10 hours/day. The follow-up period ranged from 1 to 12 years, giving a total of 2624 postoperative observations of implants, with a 3% incidence of significant skin reaction. We conclude that the surgical technique for auricular prostheses retained on osseointegrated implants is simple and associated with a low rate of postoperative and long-term complications. It offers a high degree of stability and aesthetic satisfaction. (*Otolaryngol Head Neck Surg* 1999;121:133-43)

Plastic surgical reconstruction of the defective ear is a major challenge.<sup>1</sup> In selected cases an auricular prosthesis presents an attractive and acceptable alternative. Modern techniques and plastic materials make it possible to make a prosthesis of such high quality that it is difficult to distinguish it from a normal ear.<sup>2</sup> Experience from the successful retention of dental prostheses on osseointegrated implants<sup>3</sup> has been used in developing a system for the retention of auricular prostheses. The surgical technique for the osseointegration method when used in craniofacial reconstructions is well known, and the results are well documented.<sup>4-6</sup> The biologic and physical properties of the implant system have been thoroughly investigated.<sup>7</sup>

The aim of this study was to design and present a protocol enabling systematic collection of clinical data during clinical follow-up visits to ensure the safety and quality of the system. This article describes the patients and the problems they have encountered after receiving auricular prostheses consecutively since 1979 at the Department of Otolaryngology-Head and Neck Surgery, Sahlgrenska Hospital, University of Göteborg, Sweden. The quality of the system was further investigated by a questionnaire designed to examine the way in which the prostheses affected the quality of life of the patients.

#### METHODS AND MATERIAL Protocol

The protocol presented in this study was created to enable registration (and processing) of adequate data not only for the patients but also for each individual implant (or fixture). The data registered for the patients were considered to be of importance for the outcome of the prosthesis system. The patient history includes diagnosis, general history, and when relevant, tumor data. It also includes data on smoking, medication, and diseases, along with a section on prosthesis data. The protocol data are presented as an appendix.

The protocol shows a coordinate system in which each individual fixture inserted can be given a position. The coordinates are made clockwise, with the center in the ear canal.

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### Nasal defects and osseointegrated implants: UCLA experience

Russell D. Nishimura, DDS,<sup>a</sup> Eleni Roumanas, DDS,<sup>b</sup> Peter K. Moy, DMD,<sup>c</sup> and Toshiro Sugai, DDS, PhD<sup>d</sup>  
University of California, Los Angeles, School of Dentistry, Los Angeles, Calif.

A clinical study of 33 craniofacial implants placed in 11 nasal defects was conducted over a 7-year period. Implant-retained nasal prostheses were fabricated, implant success rate was determined, and the soft tissue responses were recorded at 6-month intervals. No data were gathered on two implants because of patient death. The implant success rate was 71.4% (15/21) but varied significantly by anatomic site. The implant success rate in the glabella was 0% (0/4), whereas the success rate in the anterior nasal floor was 88.1% (15/17). All implant failures occurred within the first year of loading. A five-point scale was used to record the health of the peri-implant soft tissues, and the patients were followed up from 6 to 74 months. The unit of measure was a visit/site, and a unit was assigned for each instance an implant site was evaluated. Evaluations were conducted at 6-month intervals for a total of 76 visit/sites for the study period. The results revealed that 85.5% (65/76) of the visit/sites demonstrated an absence of inflammation; 10.5% (8/76) of the visit/sites demonstrated slight redness; 1.3% (1/76) demonstrated peri-implant red and moist tissues; 2.0% (2/76) demonstrated granulation tissue associated with the implants; and 0% (0/76) demonstrated infection of the peri-implant soft tissues. Severe soft tissue reactions around implants placed in the anterior nasal floor are rare. (J Prosthet Dent 1996;76:597-602.)

Most rhinectomy defects are acquired and are caused by resection of benign or malignant neoplasms. These defects may be surgically reconstructed or restored with a prosthesis. Surgical reconstruction of rhinectomy defects is preferable for partial nasal defects when sufficient bone and cartilaginous structures remain to support the soft tissue reconstruction. However, surgical reconstruction becomes more demanding as the defect increases in size, particularly if the patient receives cancericidal levels of radiotherapy. In these instances, prosthetic restoration is preferable because a superior esthetic result is achieved and direct visual inspection of the defect can be maintained for monitoring tumor recurrence.

In the past most nasal prostheses were retained with skin adhesives or by engaging hard or soft tissue undercuts. However, the effectiveness of adhesives was often compromised by the presence of mobile tissues in the defect (most commonly the maxillary lip), by nasal secretions, and by the warm moist air associated with respiration. The use of titanium implants for retaining na-

sal prostheses has provided significant advantages over previous methods of retention.<sup>1,2</sup>

Published implant success rates in the nasal region have been based on limited experience and have not examined the soft tissue response to the abutments. Parel and Tjellstrom<sup>3</sup> reported a 100% long-term implant success rate in Swedish patients. However, their study examined only nine implants in four patients, and none of the patients were irradiated. Granstrom et al.<sup>4</sup> continued studies on the Swedish group of patients and reported implant success rates of 50% (2/4) in the glabella and 97.5% (14/10) in the anterior nasal floor. The multicenter results compiled in the United States<sup>5</sup> and in Canada<sup>6</sup> are based on survey respondents with significantly shorter study periods than the Swedish studies. The American and Canadian implant success rates varied from 79.5% to 100%.

The purpose of this article is to report the success rates and soft tissue reactions of 23 craniofacial implants placed in the glabella and in the anterior nasal floor to retain and stabilize the nasal prostheses of 11 patients.

#### MATERIAL AND METHODS

Patients selected for the study were required to have a total rhinectomy defect and be in good physical condition. Patients with implant sites irradiated with greater than 60 Gy were excluded from the study. Additional exclusionary criteria included excessive alcohol use, drug abuse, present or latent psychotic conditions, pregnancy, and concomitant medical treatment for conditions that affected healing.

Presented at the annual meeting of the American Academy of Maxillofacial Prosthetics, Palm Springs, Calif., 1993.

This project was supported by the Albert & Elaine Berchard Foundation, Inc., and the Career Development Award, American Cancer Society.

<sup>a</sup>Assistant Professor, Section of Removable Prosthodontics.

<sup>b</sup>Adjunct Assistant Professor, Section of Removable Prosthodontics.

<sup>c</sup>Lecturer, Section of Oral and Maxillofacial Surgery.

<sup>d</sup>Visiting Associate Professor, Section of Removable Prosthodontics.

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
### Ear Epistheses as an Alternative to Autogenous Reconstruction

Philipp A. Federspil, Dr. Med.<sup>1,2</sup>

#### ABSTRACT

An ear episthesis is an artificial substitute for the auricle. The term *ear prosthesis* is used synonymously. The breakthrough came with the introduction of the modern silicones and their colorings. Although there are still indications for noninvasive methods of retention such as medical adhesives, the best and most reliable method of fixation is by bone anchorage. Long-lasting osseointegration with reaction-free skin penetration can be achieved with titanium implants. The first system used extracranially was the Hirschenmark flange fixture. Later, different solitary titanium implants were introduced, such as the ITI system. A different strategy used the titanium grids (Epitex) or plates (Epiplating) derived from osteosynthesis systems. These systems are fixed subperiosteally with several bone screws and are therefore also labeled as grouped implants. With these modern developments, secure retention can be achieved also in unfavorable anatomic situations. The grouped systems are resistant to torque with abutment insertion. The latest development is the subcutaneously implanted double magnet without skin penetration. The advantages of implant retained ear epistheses include optimal camouflage, predictable cosmetic results, fast rehabilitation, no donor site morbidity, and early detection of tumor recurrence. Depending on the clinical setting, prosthetic rehabilitation may be more than just an alternative to plastic reconstructive surgery.

**KEYWORDS:** Episthesis, auricular prosthesis, craniofacial prosthesis, titanium, osseointegrated implant, implant retained, bone anchorage, plate

An episthesis or craniofacial prosthesis is an artificial substitute for lost or removed parts of the face and adjacent structures. In the German and Scandinavian linguistic area, the word *episthesis* is preferred to stress that these types of prostheses are placed on top of some part of the face. The art of making epistheses is called *anaplastology*. Virtually all available materials e.g., porcelain, wax, rubber, papier maché have been used in the long history of anaplastology.<sup>1</sup> Obviously, the major drawbacks for rehabilitation with auricular epistheses were the use of inadequate material and the lack of

reliable methods for retention.<sup>2</sup> A major step forward was the introduction of methyl methacrylate. This material can be kept over a long lifetime and may be recolored after finishing its fabrication process. The major drawback is the stiffness of the material, which does not comply with facial mimics. However, methacrylate epistheses have retained an indication in areas where medical care cannot be given on a regular basis. The breakthrough for ear epistheses came with the introduction of the modern silicones and their colorings (Fig. 1). Silicone is flexible and keeps the body

<sup>1</sup>Department of Oto-Rhino-Laryngology, University Hospital Heidelberg, Heidelberg, Germany; <sup>2</sup>Vice president of the International Association for Surgical Prosthetics and Epistheses (IASPE, Linz, Austria).

Address for correspondence and reprint requests: Dr. med. Philipp A. Federspil, University Hospital Heidelberg, Department of Oto-Rhino-Laryngology, Im Neuenheimer Feld 405, 69120 Heidelberg,

Germany (e-mail: federspil@med.uni-heidelberg.de); Ear Reconstruction, Guest Editor, Ralf Siegent, M.D., D.D.S., Ph.D.

Facial Plast Surg 2009;25:190-203. Copyright © 2009 by Thieme Medical Publishers, Inc., 333 Seventh Avenue, New York, NY 10001, USA. Tel: +1(212) 584-4662. DOI 10.1055/s-0029-1239446. ISSN 0736-6825.

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### Treatment satisfaction with facial prostheses

Ting-Ling Chang, DDS,<sup>a</sup> Neal Garrett, PhD,<sup>b</sup> Eleni Roumanas, DDS,<sup>c</sup> and John Beumer III, DDS, MS<sup>d</sup>  
Division of Advanced Prosthodontics, Biomaterials, and Hospital Dentistry, The Weintraub Center for Reconstructive Biotechnology, UCLA School of Dentistry, Los Angeles, Calif

**Statement of problem.** Facial defects secondary to the treatment of neoplasms, congenital malformations, and trauma result in multiple functional and psychosocial difficulties. Prosthetic rehabilitation attempts to restore these facial disfigurements and may improve the level of function and self-esteem for these patients. However, a limited number of studies have evaluated the change in perceived quality of life after maxillofacial prosthetic rehabilitation.

**Purpose.** The purpose of this study was to evaluate patients' perceptions of treatment with adhesive-retained and implant-retained facial prostheses and to assess differences in overall satisfaction with these 2 types of treatments.

**Material and methods.** In this study, a questionnaire with 28 items was administered for evaluation of perceptions of appearance, comfort, fit and irritation, reliability of retention, frequency of wear, ease of placement and removal, level of self-consciousness, and value of treatment. Subjects were categorized into 2 groups: adhesive-retained group ( $n = 16$ ) and implant-retained group ( $n = 19$ ). Comparisons were made for each item in the questionnaire using Fisher exact tests ( $\alpha = .05$ ).

**Results.** The implant group reported higher positive ratings on all 28 questionnaire items when compared with the adhesive group. Statistically significant ( $P < .05$ ) differences between the implant and adhesive groups were noted for ease of placement and removal, frequency of wear at home, and quality of retention during various activities, such as home chores and when prepping or sneezing/coughing.

**Conclusion.** The implant-retained facial prosthesis offers significant enhancement over an adhesive-retained prosthesis with respect to ease of use and retention during a variety of daily activities, resulting in greater use of the prosthesis. (J Prosthet Dent 2005;94:275-80.)

### CLINICAL IMPLICATIONS

*Improvements in ease of use and retention with an implant-retained facial prosthesis appear to increase prosthesis use when compared to an adhesive-retained prosthesis. However, clinicians should evaluate patient factors, treatment costs, and the burden of additional surgery prior to determining the most appropriate prosthetic treatment for a patient.*

Facial defects secondary to the treatment of neoplasms, congenital malformations, and trauma result in multiple functional and psychosocial difficulties. Prosthetic rehabilitation to restore these facial disfigurements may improve the level of function and self-esteem for patients. However, difficulties with facial prostheses arise due to movable tissue beds, quality of prosthesis retention, and associated skin reactions to adhesives. The

use of osseointegrated implants in the craniofacial region reduces prosthesis limitations associated with medical-grade adhesives and has been proven to be a reliable treatment option with high long-term success rates for facial prostheses.<sup>1</sup> Patient acceptance of facial prostheses may be significantly enhanced due to the quality of prosthesis retention and stability afforded by craniofacial implants.

The concept of quality of life (QOL) has emerged as an organizing schema to describe and evaluate the experience of patients in clinical research. Many definitions for QOL reflect "the ability to conduct daily activities" from the patients' perspective.<sup>2</sup> There have been numerous studies reporting the QOL of head and neck cancer patients.<sup>3-10</sup> These studies indicate elevated levels of emotional distress, physical limitations, disturbed body image, and impaired relationships. Studies of the change in perceived QOL after maxillofacial prosthetic rehabilitation are limited.<sup>11-13</sup>

This investigation was conducted in a facility constructed with support from Research Facilities Improvement Program Grant No. COG-RR-14529-01 from the National Center for Research Resources, National Institutes of Health. Presented at the 6th Meeting of the International Society for Maxillofacial Rehabilitation, Maastricht, Netherlands, June, 2004.

<sup>a</sup>Clinical Associate Professor.

<sup>b</sup>Associate Professor.

<sup>c</sup>Associate Professor.

<sup>d</sup>Professor and Chair.

SEPTEMBER 2005

THE JOURNAL OF PROSTHETIC DENTISTRY 275

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